

Anti-Fibrillarin/U3 Rnp Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P22087

Item Number: CM0013389



Introduction

Rabbit polyclonal antibody targeting fibrillarin (FBL), a nucleolar methyltransferase. Validated for WB, IF/ICC, IP, and ELISA applications. Reacts with human, mouse, rat, and monkey. Ideal for ribosome biogenesis and nucleolar function research.

[Learn More](#)

Parameter	Value
Product Name	Fibrillarin/U3 RNP Rabbit pAb
Remarks/Alias	FIB; FLRN; Nop1; RNU3IP1; Fibrillarin/U3 RNP
Species	Human
GeneID (Human)	2091
GeneID	2091
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-321 of human Fibrillarin/U3 RNP (NP_001427.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	P22087
Protein Weight	34 kDa
Shipping	Ice bag