

Anti-Fibrinogen Alpha Chain (Fga) Polyclonal Antibody For Wb, Ihc-P, Elisa - P02671

Item Number: CM0018464



Introduction

Rabbit polyclonal antibody targeting human Fibrinogen alpha chain (FGA). Validated for WB, IHC-P, ELISA, cross-reacts with mouse and rat. Ideal for hemostasis and thrombosis research.

[Learn More](#)

Parameter	Value
Product Name	Fibrinogen alpha chain (FGA) Rabbit pAb
Remarks/Alias	Fib2; Fibrinogen alpha chain (FGA)
Species	Human
GeneID (Human)	2243
GeneID	2243
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 410-559 of human Fibrinogen alpha chain (FGA) (NP_000499.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02671
Protein Weight	95kDa
Shipping	Ice bag