

Anti-Fibrinogen Alpha Chain (Fga) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P02671

Item Number: CM0002407

Introduction

High-quality rabbit monoclonal antibody against human Fibrinogen alpha chain (FGA). Validated for Western blot, immunohistochemistry (paraffin), and ELISA. Cross-reacts with rat. Ideal for hemostasis, thrombosis, and blood coagulation studies.

[Learn More](#)

Parameter	Value
Product Name	Fibrinogen alpha chain (FGA) Rabbit mAb
Remarks/Alias	Fib2; Fibrinogen alpha chain (FGA)
Species	Human
Gene ID (Human)	2243
Gene ID	2243
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 700-800 of human Fibrinogen alpha chain (FGA) (P02671)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Rat
SWISS	P02671
Protein Weight	95kDa
Shipping	Ice bag