

Anti-Fibroblast Activation Protein-A (Fap) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q12884

Item Number: CM0024825



Introduction

Rabbit polyclonal antibody against human FAP (seprase, 88 kDa). Validated for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for studying ECM remodeling, wound healing, and tumor progression.

[Learn More](#)

Parameter	Value
Product Name	Fibroblast activation protein-α (FAP) Rabbit pAb
Remarks/Alias	FAPA; SIMP; DPPIV; FAPalpha; Fibroblast activation protein-α (FAP)
Species	Human
GeneID (Human)	2191
GeneID	2191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-760 of human FAP (NP_004451.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12884
Protein Weight	88kDa
Shipping	Ice bag