

Anti-Fibronectin Polyclonal Antibody For Wb, If/Icc, Elisa - P02751

Item Number: CM0014537



Introduction

High-quality rabbit polyclonal antibody against human Fibronectin for WB, IF/ICC, ELISA; cross-reacts with mouse and rat; designed for research on cell adhesion, motility, and extracellular matrix biology.

[Learn More](#)

Parameter	Value
Product Name	Fibronectin Rabbit pAb
Remarks/Alias	FN; CIG; FNZ; MSF; ED-B; FINC; GFND; LETS; GFND2; SMDCF; Fibronectin
Species	Human
GeneID (Human)	2335
GeneID	2335
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2169-2355 of human Fibronectin (NP_002017.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02751
Protein Weight	272kDa
Shipping	Ice bag