

Anti-Fibronectin Polyclonal Antibody For Wb, If/Icc, Elisa - P11276

Item Number: CM0027822



Introduction

Rabbit polyclonal antibody against fibronectin (P11276), validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Recombinant protein immunogen. For cell adhesion, fibrosis, and wound healing studies.

[Learn More](#)

Parameter	Value
Product Name	Fibronectin Rabbit pAb
Remarks/Alias	Fn; Fn-1; E330027I09
Species	Mouse
GeneID (Human)	2335
GeneID	14268
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2205-2477 of mouse Fibronectin (NP_034363.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11276
Protein Weight	273kDa
Shipping	Ice bag