

Anti-Fibronectin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P02751

Item Number: CM0018902



Introduction

Polyclonal antibody against human Fibronectin (FN1), validated for Western blot, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Suitable for cell adhesion and extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	Fibronectin Rabbit pAb
Remarks/Alias	FN; CIG; FNZ; MSF; ED-B; FINC; GFND; LETS; GFND2; SMDCF; Fibronectin
Species	Human
GeneID (Human)	2335
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02751
Protein Weight	272kDa
Shipping	Ice bag