

Anti-Fibronectin Rabbit Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P02751

Item Number: CM0006675



Introduction

High-specificity rabbit monoclonal antibody against human Fibronectin, validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Targets the extracellular matrix glycoprotein involved in cell adhesion, wound healing, and cancer. Suitable for research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	Fibronectin Rabbit mAb
Remarks/Alias	FN1; CIG; ED-B; FINC; FN; FNZ; GFND; GFND2; LETS; MSF; fibronectin; Fibronectin
Species	Human
Gene ID (Human)	2335
Gene ID	2335
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 394-614 of human Fibronectin (NP_997647.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P02751
Protein Weight	71-73kDa/221-272kDa
Shipping	Ice bag