


Anti-Fig4 Polyclonal Antibody For Wb, Elisa - Q92562

Item Number: CM0013176



Introduction

Rabbit polyclonal antibody targeting human FIG4 phosphoinositide phosphatase. Validated for Western blot and ELISA, cross-reacts with mouse. Ideal for studying PI(3,5)P2 metabolism and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	FIG4 Rabbit pAb
Remarks/Alias	BOP; YVS; BTOP; SAC3; ALS11; CMT4J; KIAA0274; dj249i4.1; FIG4
Species	Human
Gene ID (Human)	9896
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 608-907 of human FIG4 (NP_055660.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92562
Protein Weight	104kDa
Shipping	Ice bag