

Anti-Filaggrin (Flg) Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P20930

Item Number: CM0011238



Introduction

High-quality FLG Rabbit Monoclonal Antibody for IF-P, IHC-P, ELISA. Reacts specifically with human filaggrin (435kDa). Ideal for epidermal differentiation, skin barrier, and dermatological disease research. Shipped on ice bag.

[Learn More](#)

Parameter	Value
Product Name	FLG Rabbit mAb
Remarks/Alias	FLG1, ATOD2, FLG-1
Species	Human
GenelD	2312
Immunogen	Recombinant fragment
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P20930
Protein Weight	435kDa
Shipping	Ice bag