

Anti-Filaggrin (Flg) Rabbit Monoclonal Antibody For Wb, Elisa - P20930

Item Number: CM0006174



Introduction

Rabbit monoclonal antibody targeting human filaggrin (FLG), with cross-reactivity to mouse and rat. Suitable for WB and ELISA applications. Supports research on epidermal differentiation, ichthyosis, and keratinization disorders.

[Learn More](#)

Parameter	Value
Product Name	FLG Rabbit mAb
Remarks/Alias	FLG1; ATOD2; FLG-1; FLG
Species	Human
GeneID (Human)	2312
GeneID	2312
Immunogen	Recombinant protein corresponding to amino acids 1-92 of human FLG (NP_002007.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P20930
Protein Weight	435kDa
Shipping	Ice bag