


Anti-Filaggrin Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P20930

Item Number: CM0009865



Introduction

Validated rabbit monoclonal antibody targeting human filaggrin (FLG, 435 kDa) for IF-P, IHC-P, and ELISA. Filaggrin aggregates keratin filaments during epidermal differentiation, essential for skin barrier research.

[Learn More](#)

Parameter	Value
Product Name	Filaggrin Rabbit mAb
Remarks/Alias	FLG1; ATOD2; FLG-1
Species	Human
GeneID (Human)	2312
GeneID	2312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 198-288 of human FLG (NP_002007.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P20930
Protein Weight	435 kDa
Shipping	Ice bag