

Anti-Filamin A Monoclonal Antibody For Wb, Elisa - P21333

Item Number: CM0012497



Introduction

Rabbit monoclonal antibody targeting human Filamin A (FLNA), validated for WB and ELISA with mouse cross-reactivity. Ideal for studying actin cytoskeleton organization, cell migration, and signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Filamin A Rabbit mAb
Remarks/Alias	FLN; FMD; MNS; OPD; ABPX; CSBS; CVD1; FGS2; FLN1; NHBP; OPD1; OPD2; XLVD; XMVD; FLN-A; ABP-280; Filamin A
Species	Human
Gene ID (Human)	2316
Gene ID	2316
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 2548-2647 of human Filamin A (P21333).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P21333
Protein Weight	281kDa
Shipping	Ice bag