

Anti-Filamin A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P21333

Item Number: CM0022447



Introduction

Rabbit polyclonal antibody targeting human Filamin A (FLNA, ABP-280). Suitable for WB, IHC-P, IF/ICC, ELISA; cross-reacts with human, mouse, and rat. Ideal for cytoskeleton, adhesion, and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Filamin A Rabbit pAb
Remarks/Alias	FLN, FMD, MNS, OPD, ABPX, CSBS, CVD1, FGS2, FLN1, NHBP, OPD1, OPD2, XLVD, XMVD, FLN-A, ABP-280, Filamin A
Species	Human
GeneID (Human)	2316
GeneID	2316
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2348-2647 of human Filamin A (NP_001104026.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21333
Protein Weight	281kDa
Shipping	Ice bag