

Anti-Filamin-B (Flnb) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O75369

Item Number: CM0019383



Introduction

Rabbit polyclonal antibody against human FLNB (Filamin-B), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for studying actin cytoskeleton dynamics and cell migration.

[Learn More](#)

Parameter	Value
Product Name	FLNB Rabbit pAb
Remarks/Alias	AOI, FH1, SCT, TAP, LRS1, TABP, FLN-B, FLN1L, ABP-278, ABP-280, FLNB
Species	Human
Gene ID (Human)	2317
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1686-1785 of human FLNB (NP_001157789.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75369
Protein Weight	278kDa
Shipping	Ice bag