

Anti-Filamin-C Polyclonal Antibody For Wb, If/Icc, Elisa - Q14315

Item Number: CM0014138



Introduction

Rabbit polyclonal antibody against Filamin-C (FLNC) validated for Western blot, IF/ICC, and ELISA. Cross-reacts with human and mouse. Ideal for sarcomere assembly and muscle disease research.

[Learn More](#)

Parameter	Value
Product Name	FLNC Rabbit pAb
Remarks/Alias	ABPA; ABPL; FLN2; MFM5; MPD4; RCM5; CMH26; ABP-280; ABP280A; FLNC
Species	Human
GeneID (Human)	2318
GeneID	2318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2160-2340 of human FLNC (NP_001449.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14315
Protein Weight	291kDa
Shipping	Ice bag