

Anti-Filip1L Polyclonal Antibody For Wb, Elisa - Q4L180

Item Number: CM0019804

Introduction

Rabbit polyclonal antibody targeting FILIP1L (DOC1, GIP130), validated for WB and ELISA with cross-reactivity to human and mouse. Suitable for studying antiangiogenic activity and tumor vasculature.

[Learn More](#)

Parameter	Value
Product Name	FILIP1L Rabbit pAb
Remarks/Alias	DOC1; DOC-1; GIP90; GIP130; FILIP1L
Species	Human
Gene ID (Human)	11259
Gene ID	11259
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human FILIP1L (NP_878913.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q4L180
Protein Weight	130kDa
Shipping	Ice bag