

Anti-Fip1L1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q6Un15

Item Number: CM0026018



Introduction

Rabbit polyclonal antibody against human pre-mRNA 3'-end-processing factor FIP1L1. Validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Suitable for studying mRNA polyadenylation.

[Learn More](#)

Parameter	Value
Product Name	FIP1L1 Rabbit pAb
Remarks/Alias	Rhe; FIP1; hFip1; FIP1L1
Species	Human
Gene ID (Human)	81608
Gene ID	81608
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human FIP1L1 (NP_001128409.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UN15
Protein Weight	67kDa
Shipping	Ice bag