

Anti-Fkbp10 Monoclonal Antibody For Wb, If/Icc, Elisa - Q96Ay3

Item Number: CM0006511

Introduction

Rabbit monoclonal FKBP10 antibody validated for WB, IF/ICC, ELISA. Recognizes human FKBP10, a 64kDa ER-resident PPIase involved in protein folding and linked to osteogenesis imperfecta. Ideal for research and IVD development.

[Learn More](#)

Parameter	Value
Product Name	FKBP10 Rabbit mAb
Remarks/Alias	FKBP10; BRKS1; FKBP65; OI11; OI6; PPIASE; hFKBP65; FK506 binding protein 10
Species	Human
GeneID	60681
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 303-582 of human FKBP10 (NP_068758.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q96AY3
Protein Weight	39kDa/64kDa
Shipping	Ice bag