

Anti-Fkbp10 Rabbit Polyclonal Antibody For Wb - Q96Ay3

Item Number: CM0016945

Introduction

KO validated anti-FKBP10 rabbit polyclonal antibody for Western blot and ELISA, reactive with human and mouse. FKBP10 is a peptidyl-prolyl cis-trans isomerase located in the endoplasmic reticulum lumen, involved in protein folding.

[Learn More](#)

Parameter	Specification
Product Name	[KO Validated] FKBP10 Rabbit pAb
Remarks/Alias	OI6; BRKS; OI11; TLH1; BRKS1; FKBP65; PPIASE; hFKBP65; 10
Species	Human
Gene ID (Human)	60681
Gene ID	60681
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 303-582 of human FKBP10 (NP_068758.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96AY3
Protein Weight	64kDa
Shipping	Ice bag