

## Anti-Fkbp11 Polyclonal Antibody For Wb, Elisa - Q8N573

Item Number: CM0021158

### Introduction

Rabbit polyclonal antibody targeting FKBP11, suitable for Western blot and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying oxidative stress resistance, Golgi pH regulation, and neuroprotection.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | FKBP11 Rabbit pAb                                                                                                  |
| Remarks/Alias    | FKBP19; FKBP11                                                                                                     |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 51303                                                                                                              |
| GeneID           | 51303                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 28-150 of human FKBP11 (NP_057678.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q8N573                                                                                                             |
| Protein Weight   | 98kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |