

Anti-Fkbp12 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P62942

Item Number: CM0019601



Introduction

High-quality rabbit polyclonal antibody targeting human FKBP12 (FKBP1A), validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studies on TGF-beta signaling, protein folding, and calcium channel regulation.

[Learn More](#)

Parameter	Value
Product Name	FKBP12 Rabbit pAb
Remarks/Alias	FKBP1; PKC12; PKC12; FKBP12; PPIASE; FKBP-12; FKBP-1A
Species	Human
Gene ID (Human)	2280
Gene ID	2280
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 2-108 of FKBP12 (NP_000792.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62942
Protein Weight	12kDa
Shipping	Ice bag