

Anti-Fkbp1B Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P68106

Item Number: CM0021156



Introduction

FKBP1B rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat FKBP1B (12.6 kDa FK506-binding protein), involved in immunosuppression and protein folding.

[Learn More](#)

Parameter	Value
Product Name	FKBP1B Rabbit pAb
Remarks/Alias	OTK4; FKBP1L; PKBP1L; PPIase; FKBP12.6; FKBP1B
Species	Human
GeneID (human)	2281
GeneID	2281
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P68106
Protein Weight	12kDa
Shipping	Ice bag