

Anti-Fkbp2 Rabbit Polyclonal Antibody For Wb, Elisa - P26885

Item Number: CM0022559



Introduction

High-quality rabbit polyclonal antibody targeting FKBP2 (FKBP13), a peptidyl-prolyl cis-trans isomerase. Validated for Western blot and ELISA in human and mouse. Ideal for protein folding and immunophilin studies.

[Learn More](#)

Parameter	Value
Product Name	FKBP2 Rabbit pAb
Remarks/Alias	FKBP13; PPIase; FKBP-13; FKBP2
Species	Human
Gene ID (Human)	2286
Gene ID	2286
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P26885
Protein Weight	16kDa
Shipping	Ice bag