

Anti-Fkbp3 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q00688

Item Number: CM0027148



Introduction

Anti-FKBP3 rabbit polyclonal antibody validated for WB, IP, and ELISA. Reacts with human FKBP3, a nuclear immunophilin with peptidyl-prolyl isomerase activity.

[Learn More](#)

Parameter	Value
Product Name	FKBP3 Rabbit pAb
Remarks/Alias	FKBP-3; FKBP25; PPlase; FKBP-25; FKBP3
Species	Human
GeneID (Human)	2287
GeneID	2287
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q00688
Protein Weight	25kDa
Shipping	Ice bag