

## Anti-Fkbp38/Fkbp8 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q14318

Item Number: CM0010160



### Introduction

Rabbit monoclonal antibody targeting human FKBP38/FKBP8 (Q14318). Validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Useful for apoptosis, mitochondrial, and viral infection research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | FKBP38/FKBP8 Rabbit mAb                                                                                                   |
| Remarks/Alias    | FKBP38; FKBP38; FKBP38/FKBP8                                                                                              |
| Species          | Human                                                                                                                     |
| GeneID (Human)   | 23770                                                                                                                     |
| GeneID           | 23770                                                                                                                     |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FKBP38/FKBP8 (Q14318) |
| Source           | Rabbit                                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                                     |
| Application      | WB, IF/ICC, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                         |
| SWISS            | Q14318                                                                                                                    |
| Protein Weight   | 45kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |