

Anti-Fkbp38/Fkbp8 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q14318

Item Number: CM0029846



Introduction

Rabbit polyclonal antibody against human FKBP38/FKBP8, validated for WB, IF/ICC, IP, ELISA. Cross-reacts with mouse, rat. Molecular weight 45 kDa. Suitable for apoptosis and viral infection research.

[Learn More](#)

Parameter	Value
Product Name	FKBP38/FKBP8 Rabbit pAb
Remarks/Alias	FKBP38; FKBP38; FKBP38/FKBP8
Species	Human
Gene ID (Human)	23770
Gene ID	23770
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-270 of human FKBP38/FKBP8 (NP_036313.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14318
Protein Weight	45kDa
Shipping	Ice bag