


Anti-Fkbp51/Fkbp5 Polyclonal Antibody For Wb, Elisa - Q13451

Item Number: CM0027122



Introduction

Rabbit polyclonal antibody against FKBP51/FKBP5, validated for WB and ELISA in human, mouse, and rat samples. Ideal for studying steroid hormone receptor trafficking, Akt signaling, and immunophilin function.

[Learn More](#)

Parameter	Value
Product Name	FKBP51/FKBP5 Rabbit pAb
Remarks/Alias	P54; AIG6; FKBP51; FKBP54; PPIase; Ptg-10; FKBP51/FKBP5
Species	Human
GeneID (human)	2289
GeneID	2289
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 308-457 of human FKBP51/FKBP5 (NP_001139247.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13451
Protein Weight	51kDa
Shipping	Ice bag