


Anti-Fkbp51/Fkbp5 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13451

Item Number: CM0011946



Introduction

Rabbit monoclonal antibody against FKBP51/FKBP5 (Q13451), validated for WB, IHC-P, ELISA with human and rat reactivity. Target is a 51 kDa immunophilin involved in steroid receptor trafficking, Akt signaling, and NF-kB activation. Ideal for signal transduction and immunology research.

[Learn More](#)

Parameter	Value
Product Name	FKBP51/FKBP5 Rabbit mAb
Remarks/Alias	P54; AIG6; FKBP51; FKBP54; PPlase; Ptg-10; FKBP51/FKBP5
Species	Human
Gene ID (Human)	2289
Gene ID	2289
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13451
Protein Weight	51kDa
Shipping	Ice bag