

Anti-Fkbp52/Fkbp4 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q02790

Item Number: CM0027291



Introduction

Rabbit polyclonal antibody targeting human FKBP52/FKBP4, suitable for WB, IHC-P, IF/ICC, IP, ELISA. Recognizes human, mouse, and rat samples. Predicted molecular weight 52 kDa. Ideal for immunophilin and steroid receptor co-chaperone studies.

[Learn More](#)

Parameter	Value
Product Name	FKBP52/FKBP4 Rabbit pAb
Remarks/Alias	HBI; p52; Hsp56; FKBP51; FKBP52; FKBP59; PIase; FKBP52/FKBP4
Species	Human
GeneID (Human)	2288
GeneID	2288
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 220-459 of human FKBP52/FKBP52/FKBP4 (NP_002005.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02790
Protein Weight	52kDa
Shipping	Ice bag