

Anti-Fkbp52/Fkbp4 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q02790

Item Number: CM0011863



Introduction

Rabbit monoclonal antibody targeting human FKBP52/FKBP4 (SWISS: Q02790). Suitable for WB, IF/ICC, ELISA with mouse and rat cross-reactivity. Ideal for studies on steroid receptor complexes and immunophilin functions.

[Learn More](#)

Parameter	Value
Product Name	FKBP52/FKBP4 Rabbit mAb
Remarks/Alias	HBI; p52; Hsp56; FKBP51; FKBP52; FKBP59; PPIase; FKBP52/FKBP4
Species	Human
Gene ID	2288
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human FKBP4 (Q02790).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02790
Protein Weight	52kDa
Shipping	Ice bag