

Anti-Fkbp6 Polyclonal Antibody For Wb, If/Icc, Elisa - O75344

Item Number: CM0026437



Introduction

FKBP6 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets inactive peptidyl-prolyl cis-trans isomerase FKBP6 involved in spermatogenesis. Ideal for reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	FKBP6 Rabbit pAb
Remarks/Alias	FKBP36; SPGF77; FKBP6
Species	Human
GeneID (Human)	8468
GeneID	8468
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-327 of human FKBP6 (NP_003593.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75344
Protein Weight	37kDa
Shipping	Ice bag