


Anti-Fkrp Polyclonal Antibody For Wb, Elisa - Q9H9S5

Item Number: CM0013237



Introduction

FKRP rabbit polyclonal antibody validated for WB and ELISA, with cross-reactivity to rat. Targets the ribitol 5-phosphate transferase involved in dystroglycan glycosylation. Suitable for muscular dystrophy research.

[Learn More](#)

Parameter	Value
Product Name	FKRP Rabbit pAb
Remarks/Alias	FKTR; MDC1C; LGMD2I; LGMDR9; MDDGA5; MDDGB5; MDDGC5; FKRP
Species	Human
Gene ID (Human)	79147
Gene ID	79147
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-220 of human FKRP (NP_077277.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9H9S5
Protein Weight	55kDa
Shipping	Ice bag