


Anti-Flap Rabbit Monoclonal Antibody For Wb, Elisa - P20292

Item Number: CM0002168



Introduction

High-quality rabbit monoclonal antibody against human FLAP (Arachidonate 5-lipoxygenase-activating protein) for WB and ELISA applications. Validated for leukotriene biosynthesis and inflammation research. Species reactivity: human. Protein weight: 18kDa.

[Learn More](#)

Parameter	Value
Product Name	FLAP Rabbit mAb
Remarks/Alias	FLAP
Species	Human
GeneID	241 (Human)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 62-161 of human FLAP (P20292)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P20292
Protein Weight	18kDa
Shipping	Ice bag