

Anti-Flcn Polyclonal Antibody For Wb, Elisa - Q8Nfg4

Item Number: CM0025595



Introduction

Rabbit anti-FLCN polyclonal antibody validated for WB and ELISA, targeting human folliculin (64 kDa, SWISS Q8NFG4). Ideal for mTORC1 signaling, cancer, and autophagy research studies.

[Learn More](#)

Parameter	Value
Product Name	FLCN Rabbit pAb
Remarks/Alias	BHD; FLCL; DENND8B; FLCN
Species	Human
Gene ID	201163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 341-579 of human FLCN (NP_659434.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NFG4
Protein Weight	64kDa
Shipping	Ice bag