

Anti-Flg Rabbit Polyclonal Antibody For Wb, Elisa - P20930

Item Number: CM0022028



Introduction

High-quality anti-FLG rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat filaggrin, a key epidermal differentiation protein. Ideal for skin biology and keratinization research.

[Learn More](#)

Parameter	Value
Product Name	FLG Rabbit Polyclonal Antibody
Remarks / Alias	FLG1; ATOD2; FLG-1; FLG
Species	Human
Gene ID (Human)	2312
Gene ID	2312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-92 of human FLG (NP_002007.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20930
Protein Weight	435 kDa
Shipping	Ice bag