

Anti-FLI1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q01543

Item Number: CM0029524



Introduction

Rabbit polyclonal antibody against human FLI1 (EWSR2, FLI-1), validated for Western blot, immunoprecipitation, and ELISA. Recognizes FLI1 in human and mouse samples. Ideal for transcription regulation and oncogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	FLI1 Rabbit pAb
Remarks/Alias	EWSR2, FLI-1, SIC-1, BDPLT21, FLI1
Species	Human
Gene ID (Human)	2313
Gene ID	2313
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human FLI1 (NP_002008.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q01543
Protein Weight	51kDa
Shipping	Ice bag