

Anti-Flotillin-1 Monoclonal Antibody For Wb, Ihc-P, Elisa - 075955

Item Number: CM0012137



Introduction

KO-validated anti-Flotillin-1 rabbit monoclonal antibody for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets caveolar scaffolding protein FLOT1 (O75955). Ideal for membrane raft, endosome, and caveolae research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Flotillin 1 Rabbit mAb
Remarks/Alias	FLOT1; flotillin-1; Flotillin 1
Species	Human
GenelD (Human)	10211
GenelD	10211
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 155-427 of human Flotillin 1 (NP_005794.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75955
Protein Weight	47kDa
Shipping	Ice bag