

## Anti-Flotillin 1 Monoclonal Antibody For Wb, Ihc-P, Elisa - 075955

Item Number: CM0011710



### Introduction

Flotillin 1 Rabbit mAb for WB, IHC-P, ELISA. Detects human, mouse, rat FLOT1, a scaffolding protein in caveolae and membrane rafts. Supports cellular localization studies.

[Learn More](#)

| Parameter        | Value                                                                                                                         |
|------------------|-------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Flotillin 1 Rabbit mAb                                                                                                        |
| Remarks/Alias    | FLOT1; flotillin-1; Flotillin 1                                                                                               |
| Species          | Human                                                                                                                         |
| GeneID (Human)   | 10211                                                                                                                         |
| GeneID           | 10211                                                                                                                         |
| Immunogen        | Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 328-427 of human Flotillin 1 (075955).) |
| Source           | Rabbit                                                                                                                        |
| Category         | Monoclonal Antibodies                                                                                                         |
| Application      | WB, IHC-P, ELISA                                                                                                              |
| Cross Reactivity | Human, Mouse, Rat                                                                                                             |
| SWISS            | 075955                                                                                                                        |
| Protein Weight   | 47kDa                                                                                                                         |
| Shipping         | Ice bag                                                                                                                       |