

Anti-Flotillin 1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - O75955

Item Number: CM0030451



Introduction

Rabbit polyclonal antibody against Flotillin-1, validated for WB, IHC-P, IF/ICC, IP, and ELISA applications. Detects human, mouse, and rat Flotillin-1 (47 kDa). Ideal for caveolae and membrane raft research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] Flotillin 1 Rabbit pAb
Remarks/Alias	FLOT1
Species	Human
Gene ID	10211
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 155-427 of human Flotillin 1 (NP_005794.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75955
Protein Weight	47kDa
Shipping	Ice bag