

Anti-Flotillin 1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 075955

Item Number: CM0020982



Introduction

Rabbit polyclonal antibody against Flotillin 1 (FLOT1) for Western blot (WB), IHC-P, and ELISA. Validated in human, mouse, and rat. Ideal for caveolin and membrane raft research.

[Learn More](#)

Parameter	Value
Product Name	Flotillin 1 Rabbit pAb
Remarks/Alias	FLOT1; flotillin-1; Flotillin 1
Species	Human
GenelD (human)	10211
GenelD	10211
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 128-427 of human Flotillin 1 (NP_005794.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075955
Protein Weight	47kDa
Shipping	Ice bag