


Anti-Flt3 Polyclonal Antibody For Wb, Ihc-P, Elisa - P36888

Item Number: CM0018473



Introduction

Rabbit polyclonal antibody targeting human FLT3 (CD135/FLK2/STK1), validated for WB, IHC-P, ELISA in human, mouse, rat. Suitable for oncology and hematopoietic research.

[Learn More](#)

Parameter	Value
Product Name	FLT3 Rabbit pAb
Remarks/Alias	FLK2; STK1; CD135; FLK-2; FLT3
Species	Human
Gene ID	2322
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 701-800 of human FLT3 (NP_004110.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36888
Protein Weight	113kDa
Shipping	Ice bag