


Anti-Flt3/Cd135 Monoclonal Antibody For Fc, Elisa - P36888

Item Number: CM0006531



Introduction

High-quality rabbit monoclonal antibody against human Flt3/CD135, validated for flow cytometry and ELISA. Detects the FLT3 receptor tyrosine kinase, key in hematopoietic differentiation and leukemogenesis. Reacts with human samples.

[Learn More](#)

Parameter	Value
Product Name	Flt3/CD135 Rabbit mAb
Remarks/Alias	FLK2; STK1; CD135; FLK-2; Flt3/CD135
Species	Human
GeneID	2322
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-163 of human Flt3/CD135 (NP_004110.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	P36888
Protein Weight	108KDa/113kDa
Shipping	Ice bag