


Anti-Flvr2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Upi3

Item Number: CM0029402



Introduction

FLVCR2 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human and mouse FLVCR2 (CCT/SLC49A2), a choline/ethanolamine and heme transporter critical for blood-brain barrier function and mitochondrial thermogenesis.

[Learn More](#)

Parameter	Value
Product Name	FLVCR2 Rabbit pAb
Remarks/Alias	CCT; EPV; PVHH; MFSD7C; SLC49A2; C14orf58; FLVCRL14q; FLVCR2
Species	Human
GeneID (Human)	55640
GeneID	55640
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human FLVCR2 (NP_060261.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UPI3
Protein Weight	57kDa
Shipping	Ice bag