

Anti-Fmo1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q01740

Item Number: CM0024494

Introduction

High-quality rabbit polyclonal antibody against human flavin-containing monooxygenase 1 (FMO1). Validated for WB, IF-P, ELISA with human and rat cross-reactivity. Ideal for xenobiotic metabolism and taurine synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	FMO1 Rabbit pAb
Remarks/Alias	FMO1; flavin containing monooxygenase 1
Species	Human
GeneID (human)	2326
GeneID	2326
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 371-470 of human FMO1 (NP_002012.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q01740
Protein Weight	60kDa
Shipping	Ice bag