

## Anti-Fmo3 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P31513

Item Number: CM0003470

### Introduction

Rabbit monoclonal antibody against human FMO3, validated for WB, IF-P, and ELISA. Suitable for studying hepatic drug metabolism and trimethylamine N-oxide production in human, mouse, and rat samples.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | FMO3 Rabbit mAb                                                                                                      |
| Remarks/Alias    | TMAU; FMOII; dJ127D3.1; FMO3                                                                                         |
| Species          | Human                                                                                                                |
| GeneID (Human)   | 2328                                                                                                                 |
| GeneID           | 2328                                                                                                                 |
| Immunogen        | Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 300-400 of human FMO3 (P31513). |
| Source           | Rabbit                                                                                                               |
| Category         | Monoclonal Antibodies                                                                                                |
| Application      | WB, IF-P, ELISA                                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | P31513                                                                                                               |
| Protein Weight   | 60kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |