


Anti-Fmrp Monoclonal Antibody For Wb, Ihc-P, Elisa - Q06787

Item Number: CM0004303



Introduction

High-quality rabbit monoclonal antibody against human FMRP, validated for western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat. Suitable for fragile X syndrome and synaptic plasticity research.

[Learn More](#)

Parameter	Value
Product Name	FMRP Rabbit mAb
Remarks/Alias	POF, FMRP, POF1, FRAXA
Species	Human
GeneID (Human)	2332
GeneID	2332
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FMRP (Q06787).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06787
Protein Weight	71kDa
Shipping	Ice bag