


Anti-Fmrp Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q06787

Item Number: CM0027197



Introduction

Rabbit polyclonal antibody targeting FMRP (Fragile X messenger ribonucleoprotein 1). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for studies in neuronal development and synaptic plasticity.

[Learn More](#)

Parameter	Value
Product Name	FMRP Rabbit pAb
Remarks/Alias	POF; FMRP; POF1; FRAXA
Species	Human
GeneID (Human)	2332
GeneID	2332
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 500-600 of human FMRP (NP_002015.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06787
Protein Weight	71kDa
Shipping	Ice bag