

Anti-Fn1 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P02751

Item Number: CM0012185



Introduction

High-affinity rabbit monoclonal antibody targeting human fibronectin (FN1). Validated for WB, IP, IF/ICC, IHC-P, ELISA. Fibronectin is a key ECM glycoprotein regulating cell adhesion, migration, and wound healing. Ideal for diagnostic and biomedical research.

[Learn More](#)

Parameter	Value
Product Name	FN1 Rabbit mAb
Remarks/Alias	FN; CIG; FNZ; MSF; ED-B; FINC; GFND; LETS; GFND2; SMDCF
Species	Human
GeneID (Human)	2335
GeneID	2335
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2205-2477 of human FN1 (NP_997647.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P02751
Protein Weight	72 kDa/74 kDa/221-272 kDa
Shipping	Ice bag