

Anti-Fn3Krp Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ha64

Item Number: CM0015347

Introduction

Rabbit polyclonal antibody targeting FN3KRP for WB and ELISA detection in human, mouse, rat. Investigate protein deglycation pathways with a reliable ketosamine-3-kinase research reagent.

[Learn More](#)

Parameter	Value
Product Name	FN3KRP Rabbit pAb
Remarks/Alias	FN3KL; FN3KRP
Species	Human
Gene ID	79672
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-309 of human FN3KRP (NP_078895.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9HA64
Protein Weight	34kDa
Shipping	Ice bag